

## Soltis Lab CTAB DNA Extraction Protocol

(Reference: Doyle & Doyle, 1987; and Cullings 1992)

Revised November 14, 2002

1. Prepare CTAB buffer, use within 2-3 days, store capped: Add polyvinylpyrrolidone (Fisher Cat#: BP431-500) and  $\beta$ -mercaptoethanol (Fisher Cat#: BP176-100) and stir to dissolve right before starting extractions:

<u>CTAB</u>	<u>PVP</u>	<u><math>\beta</math>-merc</u>
0.5 ml	0.02g	2.5 $\mu$ l
5 ml	0.2g	25 $\mu$ l
20ml	0.8g	100 $\mu$ l

2. Weigh out **10-20 mg** of silica-dried plant tissue.
3. **Grind** tissue with blue pestles
  - a. Pestles can be reused, store in 10% bleach solution and rinse well with DI water before using them
  - b. Grinding can be done with or without the aide of liquid nitrogen and/or washed and autoclaved sand.
4. Add **500 $\mu$ l of CTAB buffer** and grind samples a bit more.
5. Incubate samples at **55°C for 1hr** to overnight.
6. Add **500 $\mu$ l of 24:1 Chloroform :Iso Amyl Alcohol** and mix well by shaking tubes.
7. Centrifuge for **5-10 minutes** at maximum speed.
  - a. Following centrifugation, you should have three layers: top: aqueous phase, middle: debris and proteins, bottom: chloroform.
  - b. Go on to the next step quickly so the phases do not remix
8. Pipette off the **aqueous phase** taking care not to suck up any of the middle or chloroform phases. Pipetting slowly helps with this.
9. Place the aqueous phase into a new labeled eppendorf tube.
10. **Estimate the volume** of the aqueous phase.
11. Add 0.08 volumes of cold **7.5 M ammonium acetate**—see attached table.
12. Add 0.54 volumes (using the combined volume of aqueous phase and added AmAc) of cold **isopropanol** (=2-propanol)—see attached table.
13. Mix well.
14. Let sit in freezer for 15 min to overnight.
  - a. Longer times will tend to yield more DNA, but also more contaminants.
15. Centrifuge for **3 min** at maximum speed.
16. Pour or pipette off the liquid, being careful not to lose the pellet with your DNA.
17. Add **700 $\mu$ l of cold 70% Ethanol** and mix
18. Centrifuge for **1 min** at maximum speed.
19. Pour or pipette off the liquid, being careful not to lose the pellet with your DNA.
20. Add **700 $\mu$ l of cold 95% Ethanol** and mix
21. Centrifuge for **1 min** at maximum speed.
22. Pour or pipette off the liquid, being careful not to lose the pellet with your DNA.
23. **Dry** the pellet:
  - a. Place samples in the speed vac for 20 min.or

- b. Invert samples on a Kim-wipe and let stand for 1 hr or until dry.
24. **Resuspend** samples with **100 $\mu$ l of TE buffer**. Allow to resuspend for 1hr at 55°C or overnight in refrigerator before running a test gel using 5 $\mu$ l of the DNA.

### **Stocks:**

**CTAB:** for 1L of CTAB buffer

100 ml of 1 M Tris, pH 8.0

280 ml of 5 M NaCl

40 ml of 0.5 M EDTA

20 g of CTAB (Cetyltrimethyl ammonium bromide, Amresco cat#:0833-1Kg)

**TE buffer:**

[Final] for 1L use:

10 mM 10 ml of 1 M Tris, pH 8.0

1 mM 2 ml of 0.5 M EDTA

**1 M Tris, pH 8.0:** for 1 L

121.1 g Tris (Fisher Cat#: BP152-5)

700 ml ddH<sub>2</sub>O

Dissolve tris and bring to 900 ml.

pH to 8.0 with concentrated HCl (will need ~50ml)

Bring to 1 L.

**0.5 M EDTA pH 8.0:** for 1 L

186.12 g of EDTA (Fisher Cat#: BP120-1)

750 ml ddH<sub>2</sub>O

Add about 20 g of NaOH pellets

Slowly add more NaOH until pH is 8.0, EDTA will not dissolve until the pH is near 8.0.

**5 M NaCl:** for 1 L

292.2 g of NaCl (Fisher Cat#: BP358-10)

700 ml ddH<sub>2</sub>O

Dissolve and bring to 1 L.

Cullings, K.W. 1992. Design and testing of a plant-specific PCR primer for ecological and evolutionary studies. *Molecular Ecology* **1**:233-240.

Doyle, J.J. and J.L. Doyle. 1987. A rapid DNA isolation procedure for small quantities of fresh leaf tissue. *Phytochemistry Bulletin* **19**:11-15.