

DNA Extraction Protocol

Modified from Edwards *et al.* (1991, NAR 19:1349)

- 1.** Place about 20 mg of plant tissue in Eppendorf tube.
- 2.** Break up tissue with pellet pestle (Kontes: Vineland, NJ).
- 3.** Add a little sterile sand and 400 μ l of extraction buffer (200 mM Tris-HCl pH 7.5, 250 mM NaCl, 25 mM EDTA, 0.5% SDS).
- 4.** Grind on ice with pellet pestle until completely ground up.
- 5.** Centrifuge in table-top microfuge at 4 °C for 2 min.
- 6.** Transfer supernatant to new Eppendorf tube (about 250 μ l).
- 7.** Add 1 volume of ice cold isopropanol and incubate on ice for 5-30 min.
- 8.** Centrifuge in microfuge at room temperature for 5 min.
- 9.** Discard supernatant and dry pellet.
- 10.** Resuspend in 100 μ l of TE, store at -20 °C.
- 11.** A 1:100 dilution of this extract should be good for PCR reactions.